Elsam III Lift Table

New Design

restraint system.

- Animals are held secure with patented sliding, no choke Elsam
- Lowest starting position of any available lift table- only 4" high. Animals can easily be walked onto the table without being lifted.
- Three 1/4" braided nylon nochoke leads supplied.
- Table is slanted 2 degrees to collect liquids, making cleaning easy.
- Wheels included for easy cleaning behind table.

Features

- Hydraulic drive for dependable, smooth, quiet operation.
- 100% Stainless Steel
- Digital weigh scale useable in all positions (1 bs. and kgs.)
- Fully automatic auto zero.
- •Bright, east to read display offers exceptional brightness at wide wiewing angles.
- Tables has raised rear lip to prevent spilling behind table.
- Three extra electrical outlets for clippers, exam light, ect.

Three-year Warranty

• Covers defects on on all parts and workmanship.

By Technidyne



TECHNIDYNE CORPORATION
2190 RT 9, UNIT D2
TOMS RIVER NJ
www.technivet.com
1-800-654-8073