

## EXPERIENCE AN ULTRA LIGHTWEIGHT WITH SMART CONTROL

# MyVet Ray W

## WEIGHT INNOVATION

(-) myvet imaging

## Wall-Mounted Dental X-Ray

- Weight innovation with CNT (Carbon Nanotube) technology
- Smart angulation for optimal result
- Secure clear images with a 0.4 mm focal spot







#### Acquire clear images with a 0.4 mm focal spot

Regardless of which intraoral sensor you use, the MyVet Ray W minimizes blur and provides additional diagnostic value with a 0.4mm focal spot.





#### LIGHTWEIGHT with CNT TECHNOLOGY



#### Smart angulations exposure compensation

Acquiring certain areas of the mouth are difficult to image, such as the rostral and caudal maxilla or the rostral mandible, due to one structure superimposed over an area of interest or difficult sensor placement. When parallel imaging techniques cannot be used, the tube head must be angled appropriately to capture these images. X-ray operators often do not change exposure settings and image quality will vary as a result of this angulation. The MyVet Ray W can optimize an image exposure setting based on the angle of the tube head.



#### Smart dial for all functions

The operating panel located on the tube head creates a simpler and faster workflow. The Smart Dial has all the setting in one dial, streamlining the workflow and decreasing preparation time. There is no need to remember complicated control panel buttons and configurations.



### **SPECIFICATIONS**

Focal Spot	0.4 mm (IEC 60336)
Tube Voltage (kV)	65 kV
Current (mA)	3.0 mA
Exposure Time Range	0.05 ~ 0.5 (0.01s increment)
Source to Skin Distance	Min. 7.9 in (200 mm)
X-Ray Field	2.4 in (60 mm) Round
Arm Length (Option of horizontal arm)	450 mm (62"), 600 mm (68"), 900 mm (80") - See note A and B for comparison
Accessories	Remote Exposure Switch
Weight	36 lbs (450 mm), 37 lbs (600 mm), 40 lbs (900mm)

\* Specifications are subject to change without prior notice.

Rayence & MyVet Inc.

Tel: 201.561.0223

81 Ruckman Rd. Unit A, Closter NJ 07624

information@myvetimaging.com www.myvetimaging.com

## DIMENSIONS (Unit: mm)



f 🛛 🔰 🗅 in 🖗

#### MyVet Ray W

myvet imaging<sup>™</sup>